

INTERNAL CONNECTION

PIESSE LINE

CATALOGUE 2023

Design & Production	05
AMS Surface Treatment	06
Package	08
Oxy Clamp	09
General Features	10
PIESSE Double Thread Micro Implants	12
PIESSE Double Thread Normo Implants	13
Prosthetic Components Healing and Transfer	14
Fixed Prosthetic Components	16
Digital Workflow	20
Cad-Cam Prosthetic Components	21
Overdenture Prosthetic Components	22
ILS Immediate Load System Prosthetic Components	24
Surgical Instruments	26
Drilling sequence	32
Guided Surgery Surgical Instruments	34

OXY
implant



OXY IMPLANT DESIGN & PRODUCTION



OXY Implant is a dental implant system totally designed, developed and manufactured in Italy.

The specific company know-how, acquired during thirty years completely spent in researches for dental implantology innovative solutions, allows to keep the **OXY Implant** line devices constantly at the quality level required by the most advanced scientific literature. All this offers the Surgeon a variety of choices that allow to easily deal with even the most complicated clinical situations.

The pursuit of the highest quality, obtained through the careful selection of suppliers and raw materials, the collaboration with research centres and the constant dialogue with the best Italian and foreign Implantologists, results in the great long-term reliability of the **OXY Implant Dental System** with the advantage of the Doctors and Patients ease.

All the implants of the **OXY Implant** line are made of Titanium Grade 4 of European origin, with characteristics of high mechanical resistance resulting from a special cold wire-drawing process.

The prosthetic components are made of Titanium Grade 5 and the surgical instruments of a special stainless steel. The machines used to produce **Oxy Implant** devices are equipped with the best numerical control technology, which allows working with tolerances of a few microns, thus ensuring the excellent overall quality of the implant-abutment system.

The production, testing and packaging phases of all the **OXY Implant** line devices are developed entirely in the company, with the consequent possibility of a direct and constant control of the whole process by a highly skilled Team of Technicians.

The primary packaging of the implants is performed in a cleanroom to avoid any contamination, in compliance with the most stringent sectoral norms.


The renovated head office allows to welcome in an appropriate and modern way all the Customers who wish to observe how the process of realisation of implants, prosthetic components and surgical instruments develops. A large training room allows to organise for Dentists and Dental Technicians many training and refresher events which also represent a useful moment of dialogue between Manufacturer and Users. Their suggestions, derived from daily practice, allow indeed to continuously improve and innovate the **OXY Implant** system.

Oxy Implant's medical devices comply with the European Regulation MDR 2017/745. The manufacturer Biomec S.r.l. has been equipped with a quality system ISO 9001 and ISO 13485 since 1998.

SURFACE TREATMENT




AMS (Advanced Micro Surface) is the surface selected for the entire range of implants by **OXY Implant**. This surface goes far beyond the results of earlier machined surfaces or surfaces treated with acid and chemical passivation. The **AMS** treatment was developed to accelerate the biological response of cell adhesion, thus enhancing the osseointegration process and the final clinical outcome. **AMS** is obtained with processes of chemical etching, decontamination and a cold Argon plasma treatment. Together they create a new surface with a high degree of cleanliness characterised by a homogeneous micro-porosity of the order of a few microns, smaller than the cell size, so as to considerably increase the adhesion of the osteoblastic cells and to favour the process of osteogenesis. All these factors are decisive in achieving the best bone-to-implant anchorage and a higher torque required for implant removal. The surface treatment of endosseous components is carried out by a specialised and certified company.



The 5000 X and 7500 X images show in detail the roughness obtained by means of the treatment: the interpeak distance is just a few microns, certainly smaller than the cell size, and accords with recent data about the effects of roughness size on the differentiation and behaviour of osteoblastic cells.

The cleanliness of the screw surface has been confirmed by an XPS analysis on the surface chemical composition. The analysed depth is about 5nm and provides a direct indication of the chemical composition of the material layers directly in contact with the bone.


According to the literature, the maximum concentration of Ti observable by means of XPS on the implant surface after treatment is 14-19%. A percentage of Titanium higher than 10% can be considered satisfactory. The detected amount of Titanium on **OXY** implants is 18%, a value close to the maximum obtainable percentage.



Cytotoxicity tests have been carried out in order to confirm the absence of toxicity induced by the processed screws. The possible presence of dead cells, multinucleated giant cells and general abnormalities in cell morphology was evaluated by comparing the data obtained at the cell layer in contact with the negative control (gold cylinder the same size as the implants) and with the positive control (gutta-percha cylinder).

After being observed under a microscope, the cells were fixed with fixative solution, coloured and photographed. The image relates to the results of tests carried out on experimental samples of the **OXY Implant** line.

All the observations reveal a situation in line with that of the negative control, that is they confirm the absence of toxic effects.



The surface treatment of the **OXY** implants guarantees:

- osteoblast adhesion to the implant surface
- no adverse effect on osteoblast adhesion and growth nor on the surface colonisation.

In conclusion:

- the surface treatment process implemented causes a morphological alteration of the processed screws
- the roughness obtained conforms to what is deemed appropriate to enhance the healing process of bone and osseointegration, at the current state of knowledge
- the implemented process and the following decontamination step allows the complete removal of residues arising from processing and the production of surfaces free of contaminants and foreign deposits
- implants processed with this treatment do not show cytotoxicity effects
- validated and checked periodically, 100% visual inspection and electron microscope analysis on some samples of each treatment batch.

PACKAGE

The implants outer package of **OXY Implant** dental system is made of cardboard and has a convenient tear-off opening system.

The packaging graphics show the indications regarding the connection type of the implant.

The external labels present the color coding and all the specific data of the implant: description, code, model number, production batch and date, expiry date, manufacturer, UDI code, CE mark, link to view the online instructions for use.

Inside the package there is the sterile tube containing the implant and the cover screw. Next to the sterile tube there are 3 stickers indicating all the specific data of the implant, which should be applied respectively on the medical record, on the patient's implant card and on any additional documentation.

An additional label is placed on the sterile tube to assure to have the main implant information available at each stage.



The prosthetic components and the surgical instruments are packaged in heat-sealed polyethylene bags, with an adhesive label showing all the data of the device: description, reference and batch numbers, color coding, manufacturer, UDI code, CE mark, link to view the online instructions for use and other specific standard symbols for each item.

The prosthetic components and surgical instruments are supplied NON-STERILE and must be sterilized before use.

The image shows four labels for the OXY Implant dental system. The top label is a white rectangular label with a QR code and UDI code (01)08054110725522 (11)230125 (10)23102722. Below the QR code, it reads 'IMP. PIESSE INT STD MICRO THREAD' and 'D. 4.50 L. 11.5' with a blue circle. The second label is a white rectangular label with a QR code and UDI code (01)08054110725522 (11)230125 (10)23102722. Below the QR code, it reads 'REF] MPSIS45115', 'LOT] 23106679', '2023-01', '2027-07', '#] 005', and 'CE 0051'. The third label is a white rectangular label with a QR code and UDI code (01)08054110725522 (11)230125 (10)23102722. Below the QR code, it reads 'IMP. PIESSE INT STD D. 4.50 L. 11.5', '#] 005', 'LOT] 23106679', 'REF] MPSIS45115', and 'CE 0051'. The fourth label is a white rectangular label with a QR code and UDI code (01)08054110725522 (11)230125 (10)23102722. Below the QR code, it reads 'IMP. PIESSE INT STD D. 4.50 L. 11.5', '#] 005', 'LOT] 23106679', 'REF] MPSIS45115', and 'CE 0051'.

OXY CLAMP PACKAGE SYSTEM

All the implants made by **OXY Implant** are contained in a plastic material tube sealed with a screw cap with safety ring, tested to guarantee a 5-year sterility.

The sterilisation is performed with BETA rays irradiation validated process, guaranteed and subjected to periodic controls.

MOUNT FREE - OXY CLAMP SYSTEM:

- implant and cover screw are placed in a plastic support with a titanium core
- it allows to see implant and cover screw before opening the sterile tube and prevents implant movements for an easier extraction

Implant extraction in 5 simple steps:



PIESSE LINE GENERAL FEATURES

The Piesse line implants, obtained from high-strength Titanium Grade 4 bars of exclusive and certified European origin, are produced by Swiss and Japanese Numerically Controlled Machines of the latest generation. Before being placed on the market, they are checked one by one by highly specialized technicians, to verify that the hexagonal connection complies with the centesimal tolerance standards established in the design phase and that the internal thread allows an easy and precise prosthetic components screwing.



Let's explore the Piesse implants more deeply:




- their morphology is conical, with different angulation degrees of the platform, central and apical areas to replicate the natural conformation of a tooth root. In particular, the apex is rounded in order to prevent the risk of anatomical structures damage during screwing. The three self-tapping longitudinal helicoidal cuttings are ideal to facilitate the insertion of the implants.
- the connection is of screwed type (M2 thread) with internal hexagon (2.41 mm key)
- the collar is machined to a height of 0.5 mm in order to minimize the probability of peri-implantitis occurrence. It has a convergent geometry also to ensure a constant aesthetic of the gingival contour, perfectly supported by the re-grown bone.
- the external thread has a double principle, in order to allow the dentist to increase the insertion speed and thus finding maximum primary stability. There are two options for the pitch: the Micro Thread of 1.2 mm and 0.6 mm intercrestal distance; and the Normo Thread, of 2.4 mm and 1.2 mm intercrestal distance.
- the external diameters are between 3.5 mm and 5 mm; the lengths instead vary between 5.5 mm and 15 mm with successive steps of 1.5 mm
- two possible connection platforms:
 - MINI** for **3.5 mm** and **4.0 mm diameter** implants with prosthetic components distinguished by the **YELLOW colour**
 - STANDARD** for implants with diameters of **4.5 mm** and **5.0 mm** with prosthetic components distinguished by the **BLUE colour**

PIESSE MICRO e NORMO THREAD

INTERNAL CONNECTION

Line: PIESSE LINE MINI
 Implant: \varnothing 3.5 - 4.0 mm
 Platform: \varnothing 3.6 mm
 ● Codice colore: YELLOW

Line: PIESSE LINE STANDARD
 Implant: \varnothing 4.5 - 5.0 mm
 Platform: \varnothing 4.0 mm
 ● Codice colore: BLUE



PIESSE DOUBLE THREAD MICRO


Ø 3.5 mm

MINI platform ø 3.6 mm
 Colour coding YELLOW
 Double Thread Micro
 Thread pitch 1.2 mm
 Intercrestal distance 0.6 mm



Ø 4.0 mm

MINI platform ø 3.6 mm
 Colour coding YELLOW
 Double Thread Micro
 Thread pitch 1.2 mm
 Intercrestal distance 0.6 mm



Ø 4.5 mm

STANDARD platform ø 4.0 mm
 Colour coding BLUE
 Double Thread Micro
 Thread pitch 1.2 mm
 Intercrestal distance 0.6 mm



Ø 5.0 mm

STANDARD platform ø 4.0 mm
 Colour coding BLUE
 Double Thread Micro
 Thread pitch 1.2 mm
 Intercrestal distance 0.6 mm



Attention: for PIESSE Micro Thread L. 5.5 mm implants is necessary to use the special green abutment screw (ref. EXDS34) supplied with the implant.

PIESSE DOUBLE THREAD NORMO

Ø 3.5 mm

MINI platform ø 3.6 mm
 Colour coding YELLOW
 Double Thread Normo
 Thread pitch 2.4 mm
 Intercrestal distance 1.2 mm



MPS35085NT



MPS35100NT



MPS35115NT



MPS35130NT



MPS35150NT

pterygoid

The Ø 4.0 and 4.5 mm pterygoid implants were created to treat upper jaw severe atrophies. The height of the machined collar is 3 mm in order to allow the tissue level placement.

Ø 4.0 mm


MINI platform ø 3.6 mm
 Colour coding YELLOW
 Double Thread Normo
 Thread pitch 2.4 mm
 Intercrestal distance 1.2 mm



MPS40070NT



MPS40085NT



MPS40100NT




MPS40115NT



MPS40130NT



MPS40150NT



MPS40170NT



MPS40200NT

Ø 4.5 mm


STANDARD platform ø 4.0 mm
 Colour coding BLUE
 Double Thread Normo
 Thread pitch 2.4 mm
 Intercrestal distance 1.2 mm



MPS45070NT



MPS45085NT



MPS45100NT



MPS45115NT



MPS45130NT



MPS45150NT



MPS45170NT



MPS45200NT

Ø 5.0 mm


STANDARD platform ø 4.0 mm
 Colour coding BLUE
 Double Thread Normo
 Thread pitch 2.4 mm
 Intercrestal distance 1.2 mm



MPS50070NT



MPS50085NT



MPS50100NT



MPS50115NT





MPS50130NT



MPS50150NT


HEALING ABUTMENT

Titanium Grade 5 | Tightening torque 10 Ncm



ANATOMICAL HEALING ABUTMENT



Titanium Grade 5 | Tightening torque 10 Ncm



Healing abutments have **different diameters and transmucosal heights**. This offers the possibility of their optimal use according to the anatomical characteristics of the gingiva. Healing abutments and transfers are also available in **anatomical** emergence profile version for optimal soft tissues management.



IMPRESSION COPING TRANSFER OPEN TRAY

Titanium Grade 5 | tightening torque 10 Ncm
impression coping transfer screw H 3 mm included ref. LTSS3402





IMPRESSION COPING TRANSFER CLOSED TRAY SNAP

Titanium Grade 5 | tightening torque 10 Ncm
impression coping transfer screw H 0 mm included ref. LTSS3400




IMPRESSION COPING TRANSFER CLOSED TRAY PULL-UP ANATOMICAL

Plastic





SNAP-ON CAP

POM



LABORATORY ANALOGUE

Stainless steel




The morphology of the **Open Tray** impression coping transfers allows perfect stability in the impression material and guarantees its absolute precision.

The particular geometry of the **Closed Tray** impression coping transfer allows maximum precision in repositioning the transfer into the silicone impression.

The analogues precisely replicate the **size** and **morphology** of the implant connection platform.


CASTABLE ABUTMENT ROTATIONAL

PMMA | Tightening torque 30 Ncm | abutment screw included




CASTABLE ABUTMENT ROTATIONAL

PMMA | Tightening torque 30 Ncm | abutment screw included




TEMPORARY ABUTMENT ROTATIONAL

Titanium Grade 5 | Tightening torque 30 Ncm | abutment screw included



TEMPORARY ABUTMENT ANTIROTATIONAL

Titanium Grade 5 | Tightening torque 30 Ncm | abutment screw included



ABUTMENT SCREW M2

Titanium Grade 5 | PACK 5 PCS



STRAIGHT ABUTMENT ANTIROTATIONAL

Titanium Grade 5 | Tightening torque 30 Ncm | abutment screw included



There are **straight and angled abutments** of traditional and anatomical shape featuring different chamfer heights and inclinations to provide a variety of solutions that deal with the most different anatomic situations. The implant-abutment connection area is realized with a 10 microns tolerance.

ANGLED ABUTMENT 15° ANTIROTATIONAL

Titanium Grade 5 | Tightening torque 30 Ncm | abutment screw included




ANGLED ABUTMENT 25° ANTIROTATIONAL

Titanium Grade 5 | Tightening torque 30 Ncm | abutment screw included



STRAIGHT ABUTMENT ANTIROTATIONAL

Titanium Grade 5 | Tightening torque 30 Ncm | abutment screw included



ANGLED ABUTMENT 15° ANTIROTATIONAL

Titanium Grade 5 | Tightening torque 30 Ncm | abutment screw included





ANGLED ABUTMENT 25° ANTIROTATIONAL

Titanium Grade 5 | Tightening torque 30 Ncm | abutment screw included



AESTHETIC ABUTMENT

Titanio Grado 5 | torque serraggio 30 Ncm | rotational castable and aesthetic abutment screws included



AESTHETIC ABUTMENT PROSTHETIC SCREW M1.8

Titanium Grade 5



The **aesthetic abutment** has a male cone that compensates for disparallelism between implants up to 55°.

It is supplied with a castable cap that allows the realization of prosthetic solutions such as: bars, Toronto Bridges, temporary abutments for immediate loading. Since the connecting screw is integrated into the abutment, the removal of the supra-structure is possible anytime without generating stress on the stability of the implant-abutment system.

SCAN BODY

opacified Titanium Grade 5 | Tightening torque 10 Ncm | screw included



SCANINT

MECHANICAL ANGLED SCREWDRIVER

Stainless steel with black diamond coating | Tightening torque 20 Ncm | torx key



SESD
short



MESD
medium

DIGITAL ANALOGUE FOR CAD/CAM

Stainless steel



ADINT

CONVERTITORE MECCANICO MANUALE

Acciaio inox temprato



MDCM8

PREMILLED ABUTMENT ANTIROTATIONAL FOR CAD-CAM

Titanium Grade 5 | abutment screw black diamond included | Arum® compatible | Ø 10 mm



INTCAD

TORX SCREW M2

Titanium Grade 5 black diamond coating



DSI3400

Implant libraries
available for:

- EXOCAD®
- 3SHAPE®
- DENTAL WINGS®

TITANIUM BASE FOR CAD-CAM

Titanium Grade 5 | Tightening torque 30 Ncm | abutment screw included




TITANIUM BASE FOR CAD-CAM

Titanium Grade 5 | Tightening torque 30 Ncm | abutment screw included





COBALT-CHROME BASE FOR CAD-CAM

Cobalt-Chrome | Tightening torque 30 Ncm | abutment screw included



BALL ABUTMENT sphere Ø 2.25 mm

Titanium Grade 5 | Tightening torque 20 Ncm | Titanium cover and retention cap included - spare part ref. 22PBEI0010T




2 pcs per package



5 pcs per package


OXY LOC ABUTMENT

Locator® compatible | Titanium Grade 5 Gold Titanium coating | Tightening torque 30 Ncm | Titanium cover, retentive cap and protective disk included (spare part ref. SPOXYLOC) | Locator® is a Zest Anchors registered trade mark

												
● mini	LOCIN2010	LOCIN2020	LOCIN2030	LOCIN2040	LOCIN2050	LOCIN2060	LOCIN2070	15LOCIN20	30LOCIN20			
												
● standard	LOCIN3410	LOCIN3420	LOCIN3430	LOCIN3440	LOCIN3450	LOCIN3460	LOCIN3470	15LOCIN34	30LOCIN34			
												
	LOCLAB laboratory analogue	LOCTR transfert	CRLOC top for angled oxy loc abutment		titanium cover for oxy loc abutment	retention cap oxy loc abutment	CN098 extra soft 0,7 kg	CN097 soft 0,9 kg	CN096 standard 1,5 kg	CN099 strong 2,8 kg	CN100 laboratory	DPR protective disk
						2 pcs per package			5 pcs per package			

STRAIGHT BASE ILS

Titanium Grade 5 | Tightening torque 30 Ncm | plastic pre-mounted handle included



ANGLED BASE 17° ILS

Titanium Grade 5 | Tightening torque 30 Ncm | stainless steel pre-mounted handle and abutment screw black diamond included



ANGLED BASE 35° ILS

Titanium Grade 5 | Tightening torque 30 Ncm | stainless steel pre-mounted handle and abutment screw black diamond included



SCREW ILS




ILASH2

ILS (Immediate Load System) includes all the components that, in a safe and reliable way, allow a fast prosthetic rehabilitation thanks to the use of All-on-Four immediate load technique.

For this purpose, **Oxy Implant** provides straight and angled bases with 17° and 35° and different transgingival heights with corresponding pre-mounted handles.


**HEALING CAP
 FOR ILS**

Titanium Grade 5 | Tightening torque 10 Ncm



**IMPRESSION COPING
 TRANSFER FOR ILS**

Titanium Grade 5 | Tightening torque 10 Ncm
 screw included



**SCAN BODY
 ILS**

Opacified Titanium Grade 5 | Tightening torque 10 Ncm
 screw included



**DIGITAL ANALOGUE
 ILS**

Stainless steel



**LABORATORY ANALOGUE
 FOR ILS**

Stainless steel



**TEMPORARY ABUTMENT
 ILS ROTATIONAL**

Titanium Grade 5 | Tightening torque 15/20 Ncm
 screw M1.4 black diamond included



**COVER
 ILS ROTATIONAL**

Tightening torque 15/20 Ncm
 screw M1.4 black diamond included



**CASTABLE
 ILS ROTATIONAL**

PMMA | Tightening torque 15/20 Ncm
 screw M1.4 black diamond included



**SCREW FOR
 PROSTHETIC ILS**

Titanium Grade 5 black diamond coating
 M1.4 | **PACK 5 PCS**



**ABUTMENT SCREW
 ILS M1.4 TORX**

Titanium Grade 5 black diamond coating | M1.4



SURGICAL INSTRUMENTS

The drills for the **PIESSE Line** implants are designed and manufactured with the aim of ensuring their constant efficiency and maximum longevity. These prerogatives are achieved first of all thanks to the use of top quality materials such as the special tempered stainless steel used for their production. The latest generation finishing processes complete the production cycle.



All drills surface is coated with **Black Diamond**, a perfectly biocompatible material which increases their surface hardness and decreases the generated frictional forces. This process therefore greatly increases the cutting power and consequently reduces the risk of bone overheating during surgery.

Their morphology is cylindrical with a step of lower diameter to make a surgical alveolus that faithfully follows the shape of **Piesse** implant, more conical in the apical area.

The drills are sequential with a step of 0.3 mm between the successive diameters (2.9 - 3.2 - 3.5 - 3.8 and 4.1 mm) to make a surgical alveolus that corresponds perfectly to the implant core or to underprepare or over-prepare it according to the characteristics of the bone density.

In order to ease the Surgeon's work, the drills are graduated with laser markings to indicate the length of the inserted implant. For the same reason they can also be equipped with a Titanium depth stopper.

As well the diameter of the drill is clearly highlighted thanks to a specific laser marking.

The colour of the instrument stem allows to easily identify the slot of the surgical kit in which to store it: it is indeed characterized by the same coding.

It is important to remember that the drill tip increases the length of the instrument by 0.5 mm. Considering this, it will always be possible to avoid damage to the anatomical structures especially during use in areas close to the maxillary sinus or to the mandibular canal.



**PILOT
DRILL**

∅ 1.8 mm

Hardened stainless steel
black diamond coating

max 900 Rpm



SID01L

**CALIBRATED
STOP**

Titanium Grade 5



Length	Ref.
5.5 mm	EST4055T
7 mm	EST4070T
8.5 mm	EST4085T
10 mm	EST4100T
11.5 mm	EST4115T
13 mm	EST4130T
15 mm	EST4150T

**GRADUATED
DRILL**

∅ 2.3 mm

Hardened stainless steel
black diamond coating

colour coding: white

max 900 Rpm



IGSD230

**GRADUATED
DRILL**

∅ 2.6 mm

Hardened stainless steel
black diamond coating

colour coding: purple

max 900 Rpm



SD2620

**GRADUATED
DRILL**

∅ 2.9 mm

Hardened stainless steel
black diamond coating

colour coding: yellow

max 900 Rpm



SD2923

**GRADUATED
DRILL**

∅ 3.2 mm

Hardened stainless steel
black diamond coating

colour coding: red

max 800 Rpm



SD3226

**GRADUATED
DRILL**

∅ 3.5 mm

Hardened stainless steel
black diamond coating

colour coding: green

max 800 Rpm



SD3529

**GRADUATED
DRILL**

∅ 3.8 mm

Hardened stainless steel
black diamond coating

colour coding: blue

max 700 Rpm



SD3832

**GRADUATED
DRILL**

∅ 4.1 mm

Hardened stainless steel
black diamond coating

colour coding: brown

max 700 Rpm



SD4135

LONG GRADUATED DRILL

PTERYGOID IMPLANT DRILLS

Hardened stainless steel
black diamond coating

- Ø 2.3 max 900 Rpm
- Ø 2.6 max 900 Rpm
- Ø 2.9 max 900 Rpm
- Ø 3.2 max 800 Rpm



LD2320
Ø 2.3 mm



LD2620
Ø 2.6 mm



LD2923
Ø 2.9 mm



LD3226
Ø 3.2 mm

TISSUE PUNCH

Hardened stainless steel
black diamond coating

max 100 Rpm



MPU34

COUNTERSINK DRILL

Hardened stainless steel
black diamond coating

max 300 Rpm



STDPSILS

PARALLEL INDICATOR

Titanium Grade 5



PIN2329
Ø2.3-2.9 mm



PIN2317
17°



PIN2330
30°

DRILL EXTENSION

Hardened stainless steel



DEXT

BONE TAP

Hardened stainless steel
black diamond coating

max 40 Rpm



double thread **micro**

Diameter	Ref.
3.5 mm	TDLT350
4.0 mm	TDLT400
4.5 mm	TDLT450
5.0 mm	TDLT500



double thread **normo**

Diameter	Ref.
3.5 mm	TNT350
4.0 mm	TNT400
4.5 mm	TNT450
5.0 mm	TNT500

BONE TAP DRIVER

Hardened stainless steel



XCMRC8
extrashort



CMRC8
short



LMRC8
medium

MECHANICAL INSERTER


Hardened stainless steel
black diamond coating
Tightening torque 50 Ncm



SMEIT short
MMEIT medium
LMEIT long

MANUAL DRIVER INSERTER

disc: INT
Hardened stainless steel
Ø 8 mm



SMI8 short
MMI8 medium
LMI8 long
EMI8 extralong

MECHANICAL SCREWDRIVER HEX 1.25


Hardened stainless steel
black diamond coating
hexagonal key 1.25 mm



MESSD short
MEMSD medium
MELSD long
MEELSD extralong

MANUAL SCREWDRIVER HEX 1.25

disc: 1.25
Hardened stainless steel
gold titanium coating
hexagonal key 1.25 mm



SMSD8 short
MMSD8 medium
LMSD8 long
ELMSD8 extralong

SCREWDRIVER FOR ILS STRAIGHT BASE AND AESTHETIC ABUTMENT

rondella: ILS
Hardened stainless steel
gold titanium coating



MADD8
MIAL

SCREWDRIVER FOR BALL ABUTMENTS AND MICRO FIX SPHERE IMPLANTS

rondella: SPHERE
Acciaio inox temprato



MBDD8
MAXPS

SURGICAL DIRECTION GUIDE ILS

Anodised Titanium



ILSGUIDE

CYLINDRICAL BONE MILL

Ø 5.0 mm
Hardened stainless steel
black diamond coating
max 300 Rpm



BMCIL



GBMSINT
bone mill screw

CONICAL BONE MILL

Ø 6.5 mm
Hardened stainless steel
black diamond coating
max 300 Rpm



BMCON



GBMSINT
bone mill screw

DEPTH GAUGE GRADUATED

Stainless steel
with depth indicator and ruler



DIPSK1

TORQUE RATCHET


Hardened stainless steel
driver D8
torque 10÷70 Ncm



TR8

STRAIGHT KEY


Hardened stainless steel




DEK

PIESSE LINE SEQUENZA FRESE

for \varnothing 3.5 mm implants




for \varnothing 4.0 mm implants




PIESSE LINE DRILLING SEQUENCE

for \varnothing 4.5 mm implants



for \varnothing 5.0 mm implants



OXY GUIDED SURGERY KIT



OXY GUIDED SURGERY **SURGICAL INSTRUMENTS**

DRILL FOR FIXATION PIN

Hardened stainless steel
black diamond coating
max 1000 Rpm



DPFGS

FIXATION PIN FOR SURGICAL GUIDES

Stainless steel



PINGS

SCREWABLE



PINGF

OGS REDUCTION RING

Titanium Grade 5



Height	Colour	Ref.
2 mm	fuchsia	RR5020
1.5 mm	blue	RR5015
1 mm	yellow	RR5010

OGS BONE LEVELLER DRILL

Hardened stainless steel
black diamond coating

max 300 Rpm



OGSBL

OGS START DRILL

Hardened stainless steel
black diamond coating
max 600 Rpm



CSGS

OGS DRILL

Ø 2.3 mm

Hardened stainless steel
black diamond coating

max 900 Rpm



Length	Ref.
L. 6.0 mm	SD23060GS
L. 8.5 mm	SD23085GS
L. 10 mm	SD23100GS
L. 11.5 mm	SD23115GS
L. 13 mm	SD23130GS
L. 15 mm	SD23150GS

OGS DRILL

Ø 2.9 mm

Hardened stainless steel
black diamond coating

max 900 Rpm



Length	Ref.
L. 8.5 mm	SD29085GS
L. 10 mm	SD29100GS
L. 11.5 mm	SD29115GS
L. 13 mm	SD29130GS
L. 15 mm	SD29150GS

OGS DRILL

Ø 3.2 mm

Hardened stainless steel
black diamond coating

max 800 Rpm



Length	Ref.
L. 8.5 mm	SD32085GS
L. 10 mm	SD32100GS
L. 11.5 mm	SD32115GS
L. 13 mm	SD32130GS
L. 15 mm	SD32150GS

OGS DRILL

Ø 3.5 mm

Hardened stainless steel
black diamond coating

max 800 Rpm



Length	Ref.
L. 8.5 mm	SD35085GS
L. 10 mm	SD35100GS
L. 11.5 mm	SD35115GS
L. 13 mm	SD35130GS
L. 15 mm	SD35150GS

OGS DRILL

Ø 3.8 mm

Hardened stainless steel
black diamond coating

max 700 Rpm



Length	Ref.
L. 8.5 mm	SD38085GS
L. 10 mm	SD38100GS
L. 11.5 mm	SD38115GS
L. 13 mm	SD38130GS
L. 15 mm	SD38150GS

OXY GUIDED SURGERY **SURGICAL INSTRUMENTS**

**MOUNTER
OGS**

Hardened stainless steel



GMIS

**MONTER
REMOVAL
DRIVER OGS**

Hardened stainless steel



GMSE

**DRIVER
CONVERTER
4X4 SQUARE**

Hardened stainless steel

manual

mechanical



SMRC8



MAEX

**CYLINDRICAL
BONE MILL**

Hardened stainless steel
black diamond coating

max 300 Rpm



BMCIL



GBMSINT
bone mill screw

**CONICAL
BONE MILL**

Hardened stainless steel
black diamond coating

max 300 Rpm



BMCON



GBMSINT
bone mill screw

**SCREWDRIVER
FOR SCREWS**

Hardened stainless steel
black diamond coating

max 300 Rpm

manual

mechanical



MMSD8



MEMSD

**SCREWDRIVER
FOR ILS
STRAIGHT BASE
AND AESTHETIC
ABUTMENT**

rondella: ILS

Hardened stainless steel
gold titanium coating

manual

mechanical



MADD8



MIAL

OXY GUIDED SURGERY **SURGICAL INSTRUMENTS**

**SLEEVE
FOR
TEMPLATE**

Stainless steel



CGS

**FIXING PIN
SLEEVE**

Stainless steel



CPINGS

**BITE
GAUGE**

Stainless steel



BITEGAUGEPS

**TORQUE
RATCHET**

Stainless steel

driver D8

torque 10÷70 Ncm



TR8




By Biomec S.r.l.

Via Nazionale Nord, 21/A - 23823 Colico (LC) - Italy

Tel. +39 0341 930166 - Fax +39 0341 930201

www.oxyimplant.com - info@oxyimplant.com

DESIGN & PRODUCTION
 100% OXY IMPLANT
MADE IN ITALY